

Бюретки

Описание

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: nxa@nt-rt.ru || сайт: <https://normax.nt-rt.ru/>



Burette w/automatic zero point two stopcock glass key Cl.AS 0,1 25 ml

Ref: 4.3431833

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... Glass



Burette w/automatic zero point two stopcock PTFE key Cl.AS 0,1 50 ml

Ref: 4.3451836

Capacity..... 50 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... PTFE





Automatic burette acc. Dr. Schilling w/ 500 ml flask class AS divis. 0,02 5 ml
Ref: 3.3086922

Capacity..... 5 ml

Certification..... ISO 385

Class Class AS

Material..... Borosilicate



Burette w/automatic zero point amber glass two stopc. PTFE key Cl.AS 0,1 25 ml

Ref: 4.3461833





Burette w/automatic zero point two stopcock PTFE key Cl.AS 0,02 10 ml

Ref: 4.3451822

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... PTFE





Burette w/automatic zero point two stopcock glass key Cl.AS 0,05 25 ml

Ref: 4.3431834

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... Glass



Burette w/automatic zero point ambar glass two stopc. glass key Cl.AS 0,05 10 ml

Ref: 4.3441827

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Colour..... Amber

Print enamel..... White

Material..... Borosilicate

Tolerance..... $\pm 0,02$ ml

Stopcock..... Glass



Automatic burette acc. Dr. Schilling w/ 1000 ml flask class AS divis. 0,1 25 ml

Ref: 3.3086933

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Material..... Borosilicate



Burette w/automatic zero point two stopcock PTFE key Cl.AS 0,05 25 ml

Ref: 4.3451834

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... PTFE





Automatic burette acc. Dr. Schilling w/ 500 ml flask class AS divis. 0,05 10 ml

Ref: 3.3086927

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Material..... Borosilicate



Burette w/automatic zero point two stopcock PTFE key Cl.AS 0,1 25 ml

Ref: 4.3451833

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... PTFE





Burette w/automatic zero point ambar glass two stopc. glass key Cl.AS 0,1 50 ml

Ref: 4.3441836

Capacity..... 50 ml

Certification..... ISO 385

Class Class AS

Colour..... Amber

Print enamel..... White

Material..... Borosilicate

Tolerance..... $\pm 0,05$ ml

Stopcock..... Glass



Burette w/automatic zero point ambar glass two stopc. glass key Cl.AS 0,1 25 ml

Ref: 4.3441833

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Colour..... Amber

Print enamel..... White

Material..... Borosilicate

Tolerance..... $\pm 0,05$ ml

Stopcock..... Glass



Burette w/automatic zero point two stopcock PTFE key Cl.AS 0,05 10 ml

Ref: 4.3451827

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... PTFE





Automatic burette acc. Dr. Schilling w/ 1000 ml flask class AS divis. 0,1 50 ml
Ref: 3.3086936

Capacity..... 50 ml

Certification..... ISO 385

Class Class AS

Material..... Borosilicate



Burette w/automatic zero point two stopcock glass key Cl.AS 0,05 10 ml

Ref: 4.3431827

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... Glass



Burette w/automatic zero point two stopcock glass key Cl.AS 0,1 50 ml

Ref: 4.3431836

Capacity..... 50 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... Glass



Burette w/automatic zero point amber glass two stopc. glass key Cl.AS 0,05 25 ml

Ref: 4.3441834

Capacity..... 25 ml

Class Class AS

Colour..... Amber

Print enamel..... White

Material..... Borosilicate

Tolerance $\pm 0,03$ ml

Stopcock..... Glass



Burette w/automatic zero point two stopcock glass key Cl.AS 0,02 10 ml

Ref: 4.3431822

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Colour..... Clear

Material..... Borosilicate

Stopcock..... Glass



Burette straight stopcock. PTFE key Cl. AS 0,2 100 ml

Ref: 4.5433039

Capacity..... 100 ml

Certification..... ISO 385

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,1$ ml

Stopcock..... PTFE



Burette straight stopcock. PTFE key Cl. AS 0,02 5 ml

Ref: 4.5433022

Capacity..... 5 ml

Certification..... ISO 385

Class Class AS

Material..... Borosilicate

Tolerance..... $\pm 0,01$

Stopcock..... PTFE



Burette straight stopcock. PTFE key Cl. AS 0,05 25 ml

Ref: 4.5433034

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,03$ ml



Burette straight stopcock. PTFE key Cl. AS 0,02 10 ml

Ref: 4.5433027

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,02$ ml



Burette straight stopcock. PTFE key Cl. AS 0,1 25 ml

Ref: 4.5433033

Capacity..... 25 ml

Certification..... ISO 385

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,05$ ml

Stopcock..... PTFE



Burette straight stopcock. PTFE key Cl. AS 0,1 50 ml

Ref: 4.5433036

Capacity..... 50 ml

Certification..... ISO 385

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,05$ ml

Stopcock..... PTFE



Burette w/automatic zero point amber glass two stopc. PTFE key Cl.AS 0,1 50 ml

Ref: 4.3461836



Burette straight stopcock. PTFE key with flared top 0,01 Cl. AS 2 ml

Ref: 4.5473016



Burette straight stopcock, PTFE key with flared top 0,2 Cl. AS 100 ml

Ref: 4.5473039



Burette straight stopcock, PTFE key flared top w/Schellbach 0,01 Cl. AS 2 ml

Ref: 4.54730161



Burette straight stopcock, PTFE key flared top w/Schellbach 0,02 Cl. AS 10 ml

Ref: 4.54730271



Burette straight stopcock, glass key Cl. AS 0,01 2 ml

Ref: 4.5432816



Burette straight stopcock, PTFE key Cl. AS 0,05 10 ml

Ref: 4.5433029

Capacity..... 10 ml

Certification..... ISO 385

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,02$ ml

Stopcock..... PTFE



Burette straight stopcock, glass key with Schellbach Cl. AS 0,01 2 ml

Ref: 4.54328161



Burette straight stopcock, PTFE key with flared top 0,1 Cl. AS 50 ml

Ref: 4.5473036



Burette straight stopcock, glass key with flared top 0,02 Cl. AS 10 ml

Ref: 4.5472827



Burette straight stopcock, glass key flared top with Schellbach 0,1 Cl. AS 50 ml

Ref: 4.54728361



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,05 10 ml

Ref: 4.54330291



Burette straight stopcock glass key Cl. AS 0,02 5 ml

Ref: 4.5432822



Burette straight stopcock, PTFE key flared top w/Schellbach 0,02 CL. AS 5 ml

Ref: 4.54730221



Burette straight stopcock. PTFE key with flared top 0,1 AS 25 ml

Ref: 4.5473033



Burette straight stopcock. PTFE key with flared top 0,02 Cl. AS 5 ml

Ref: 4.5473022



Automatic burette with two-way stopcock, PTFE key Cl. AS 0,1 50 ml

Ref: 4.3450236



Burette straight stopcock, glass key with Schellbach Cl. AS 0,05 10 ml

Ref: 4.54328291



Burette straight stopcock, PTFE key flared top w/Schellbach 0,1 Cl. AS 50 ml

Ref: 4.54730361



Burette straight stopcock, glass key with Schellbach Cl. AS 0,2 100 ml

Ref: 4.54328391



Burette straight stopcock, glass key flared top with Schellbach 0,05 Cl. AS 25 m

Ref: 4.54728341



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,2 100 ml

Ref: 4.54330391



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,02 5 ml

Ref: 4.54330221



Burette straight stopcock, glass key flared top with Schellbach 0,1 Cl. AS 25 ml

Ref: 4.54728331



Burette straight stopcock, PTFE key flared top 0,05 Cl. AS 10 ml

Ref: 4.5473029



Burette straight stopcock, PTFE key flared top w/Schellbach 0,05 Cl. AS 25 ml

Ref: 4.54730341



Burette straight stopcock glass key with flared top 0,05 Cl. AS 25 ml

Ref: 4.5472834



Burette straight stopcock, glass key with flared top 0,05 Cl. AS 10 ml

Ref: 4.5472829



Burette straight stopcock. PTFE key with flared top 0,02 Cl. AS 10 ml

Ref: 4.5473027



Burette straight stopcock, PTFE key flared top w/Schellbach 0,05 Cl. AS 10 ml

Ref: 4.54730291



Burette straight stopcock, glass key with flared top 0,01 Cl. AS 2 ml

Ref: 4.5472816



Burette straight stopcock glass key Cl. AS 0,1 50 ml

Ref: 4.5432836



Burette straight stopcock glass key Cl. AS 0,02 10 ml

Ref: 4.5432827



Burette straight stopcock, glass key Cl. AS 0,05 10 ml

Ref: 4.5432829



Burette straight stopcock. PTFE key Cl. AS 0,01 2 ml

Ref: 4.5433016

Capacity..... 2 ml

Class Class AS

Print enamel..... Blue

Material..... Borosilicate

Tolerance..... $\pm 0,01$

Stopcock..... PTFE



Burette straight stopcock, glass key with flared top 0,1 Cl. AS 25 ml

Ref: 4.5472833



Burette straight stopcock, glass key with Schellbach CL. AS 0,02 10 ml

Ref: 4.54328271



Burette straight stopcock, glass key with flared top 0,02 Cl. AS 5 ml

Ref: 4.5472822



Burette straight stopcock, glass key with Schellbach Cl. AS 0,1 50 ml

Ref: 4.54328361



Burette straight stopcock glass key with flared top 0,2 Cl. AS 100 ml

Ref: 4.5472839



Burette straight stopcock glass key Cl. AS 0,2 100 ml

Ref: 4.5432839



Burette straight stopcock, glass key with Schellbach Cl. AS 0,02 5 ml

Ref: 4.54328221



Burette straight stopcock glass key Cl. AS 0,05 25 ml

Ref: 4.5432834



Burette straight stopcock, glass key with Schellbach 0,1 Cl. AS 25 ml

Ref: 4.54328331



Burette w/automatic zero point amber glass two stopc. PTFE key Cl.AS 0,05 10 ml

Ref: 4.3461827



Burette straight stopcock, PTFE key with flared top 0,05 Cl. AS 25 ml

Ref: 4.5473034



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,1 50 ml

Ref: 4.54330361



Burette straight stopcock glass key Cl. AS 0,1 25 ml

Ref: 4.5432833



Burette straight stopcock glass key with flared top 0,1 Cl. AS 50 ml

Ref: 4.5472836



Burette straight stopcock, PTFE key flared top with Schellbach 0,1 AS 25 ml

Ref: 4.54730331



Burette straight stopcock, PTFE key flared top w/Schellbach 0,2 Cl. AS 100 ml

Ref: 4.54730391



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,02 10 ml

Ref: 4.54330271



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,05 25 ml

Ref: 4.54330341



Burette straight stopcock, PTFE key with Schellbach Cl. AS 0,1 25 ml

Ref: 4.54330331



Burette straight stopcock, glass key with Schellbach Cl. AS 0,05 25 ml

Ref: 4.54328341



По вопросам продаж и поддержки обращайтесь:				
Алматы (727)345-47-04 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волгоград (844)278-03-48 Вологда (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89	Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81	Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73 Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47	Ростов-на-Дону (863)308-18-15 Рязань (4912)46-61-64 Самара (846)206-03-16 Санкт-Петербург (812)309-46-40 Саратов (845)249-38-78 Севастополь (8692)22-31-93 Саранск (8342)22-96-24 Симферополь (3652)67-13-56 Смоленск (4812)29-41-54 Сочи (862)225-72-31 Ставрополь (8652)20-65-13 Сургут (3462)77-98-35 Сыктывкар (8212)25-95-17 Тамбов (4752)50-40-97 Тверь (4822)63-31-35	Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удэ (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93
Россия +7(495)268-04-70	Казахстан +7(727)345-47-04	Беларусь +(375)257-127-884	Узбекистан +998(71)205-18-59	Киргизия +996(312)96-26-47

эл.почта: nxa@nt-rt.ru || сайт: <https://normax.nt-rt.ru/>